
## If someone collapses, follow these instructions



**The Community HeartBeat Trust** 



## ADDITIONAL COMPRESSION/BREATHING INFORMATION

- If you are unable or unwilling to give rescue breaths just do the compressions continuously.
- Give 30 compressions at a rate of 100-120 per minute and then give 2 breaths.
- Repeat 30 compressions and 2 breaths sequence.
- Continue until help arrives, you are too tired to continue, or the casualty starts breathing normally again.

FOR AN ADULT USE 2 HANDS & COMPRESS 5-6 cm

FOR A CHILD USE 1 HAND & COMPRESS 1/3 CHEST DEPTH
FOR AN INFANT USE 2 FINGERS & COMPRESS 1/3 CHEST DEPTH

TO THE COUNTY OF THE COUNTY OF

Augmented reality – if you would like to see some videos on how to do chest compressions and how to use a defibrillator, use your smart phone to scan these QR codes.

Community Heartbeat Trust



CPR



Defibrillator



The Community Heartbeat Trust – the UK's community defibrillator charity – www.communityheartbeat.org.uk